
The Monitoring Board for the Calibration System of the Muon g-2 Experiment

Octavio Escalante New Perspectives 2016 Fermilab | 14 June 2016







Overview



- Experiment Background
- Laser Calibration System
- Monitoring Board
- Summary


Experiment Background: "Precision" frontier







Does particle behavior match Standard Model predictions?

- Mass?
- Production rates?
- Decay rates?
- Interactions with other particles or fields (e.g. magnetic moment)




Experiment Background: "Precision" frontier






Does particle behavior match Standard Model predictions?

- Mass?
- Production rates?
- Decay rates?
- Interactions with other particles or fields (e.g. magnetic moment)



Experiment Background: Magnetic moment





QED


$$\vec{\mu} = g \; \frac{q}{2m} \vec{S}$$

Dirac theory: a charged, spin ½ elementary point particle has

$$g \equiv 2$$








$$a_{\mu}^{SM} = (g_{\mu}^{SM}-2)/2 = a_{\mu}^{QED} + a_{\mu}^{EW} + a_{\mu}^{QCD}$$


Experiment Background: Testing the anomalous magnetic moment




$$a_{\mu}^{\text{Expt.}}$$
 - a_{μ}^{SM} = $(260 \pm 78) \times 10^{-11}$ (3.3 σ)
$$a_{\mu}^{\text{Expt.}} = a_{\mu}^{\text{SM}} + a_{\mu}^{\text{New Physics}}$$


- E-821 at BNL
 - Latest measurement of the anomalous magnetic moment of a muon had a 3.3σ discrepancy from SM
 - Uncertainty mainly due to HVP and HLBL terms in QCD prediction
- E-989 at Fermilab
 - More than 21 times the amount of statistics than predecessor E-821
 - δa_μ^{exp} = .54 ppm to .14 ppm improvement
 - Reduced pion contamination, segmented detectors and an improved storage ring kicker

Experiment Background: Muons in a storage ring




- 1. Start with polarized muon beam (from pion decay)
- 2. Cyclotron frequency: $\omega_c = \frac{e}{m \gamma} B$
- 3. Spin precession frequency: $\omega_{\rm S} = \frac{e}{m \; \gamma} \; B \; (1 + \gamma \; a_{\mu})$






$$\omega_a = \omega_S - \omega_C = e/m a_{\mu} B$$




Jegerlehner & Nyffeler, Phys. Rept. 477 (2009) 1-110, arXiv:0902.3360v1

Experiment Background: Muons in a storage ring



$$N(t) = N_0 e^{-t/\tau} \left(1 + A\cos(\omega_a t + \phi) \right)$$

au is the muon lifetime, the magnitude of the amplitude A is determined by the energy cut, and the phase ϕ depends on the initial polarization of the muon ensemble.




- Self-analyzing decay:
 - Higher energy positrons emitted preferentially in direction of muon spin
- Spectrum distortions
 - Pileup, gain stability
 - Beam effects, losses


Time spectrum of decay positrons above 1.8 GeV/c^2 . Modulation is at frequency ω_a which is proportional to a_μ (Courtesy of the E821 collaboration)


Experiment Background: Instrumentation Upgrades
















Experiment Background: Instrumentation Upgrades










Segmented compact PbF2 Cherenkov







Laser Calibration System: Variable gain





Gain fluctuations


– Dangerous "short term" (within 700 µs fill) fluctuations alter energy reconstruction and raise systematic error on ω_a to intolerable levels (>.01 ppm)

Calibration approach


- Monitor gain (G_{cal}) of all calorimeter elements by exciting and monitoring the entire system periodically with a <u>common</u> light source (in principle, can ignore light source fluctuations)
- Statistical fluctuations in recorded data must be smaller than variations in G_{cal} (<.1% per hour)








Laser Calibration System: SM detectors




Each Source Monitor detects light using 3 independent detectors, 2 PIN diodes and 1 PMT to observe for eventual beam pointing effects.

- PIN photodiodes (S1722-02)
 - amplified by custom frontend electronics
 - High photoefficiency (> 70%)
 - Fast (rise time ~ ns) can be shaped according to necessity
 - Stable inherent gain (rise time and amplitude depend on bias but charge is nearly independent of bias voltage)
 - Overall detector gain stability depends essentially on frontend electronics which are designed to optmize stability


Laser Calibration System: SM detectors





- PMT (H5783-04)
 - amplified by custom frontend electronics
 - receives light pulses transmitted from the mixing chamber to the photocathode by wavelength-shifting fibers
 - light pulses which are emitted by a weak (~ 5-10 Hz) Am source (deposited on an Nal crystal, enclosed in an aluminium cylinder with a quartz window) are situated close to the PMT photocathode
 - the signal from the Am source serves as an absolute reference to aid with the relatively poor stability (e.g. strong dependence on HV) of PMTs



UniversiTàdegli STudi di Napoli Federico II








UniversiTàdegli STUDI di Napoli Federico II

Input pulse: 180 pJ


Run time: 16 h


Pulse rate: < 1 Hz,

WFD: DRS4




Chanel output (V) vs Time (ns)






Closing remarks




- Storage ring is cold and powered
 - Installation in 2nd half of 2016
 - Muons in 2017
- Laser calibration system will aid in obtaining 0.01 ppm statistical uncertainty
 - Monitoring equipment is ready for installation
 - Temperature and electronic baseline fluctuations must be well understood (frontend electronics)







Backup




Experiment Background: Measuring aµ






How it works


- 3.1 GeV/c muons produced from pion decay channel are fed into 14m diameter, 1.45T magnet
- Decaying muons emit electrons which curl towards ring interior and are collected by calorimeters







PHOTON IS OUR BUSINESS



Si PIN photodiode

S3590-08/-09/-18/-19

Large active area Si PIN photodiode

➡ Structure / Absolute maximum ratings


Type No.			Daulatian	Absolute maximum ratings					
	Window material	Active area	Depletion layer thickness (mm)	Reverse voltage VR max	Power dissipation P (mW)	Operating temperature Topr (°C)	Storage temperature Tstg (°C)		
S3590-08	Epoxy resin		0.3	100	100	-20 to +60	20 +- + 00		
S3590-09	Unsealed	10 × 10							
S3590-18	Epoxy resin				100		-20 to +80		
S3590-19	Unsealed	Ī							

Note: Exceeding the absolute maximum ratings even momentarily may cause a drop in product quality. Always be sure to use the product within the absolute maximum ratings.

₽ Electrical and optical characteristics (Typ. Ta=25 °C, unless otherwise noted)

Type No.	Spectral response range λ	Peak sensitivity wavelength λp	Photo sensitivity S				Short circuit current	uit currer ID ent VR=70		nt lemp. coefficient of ID	Cut-off Frequency fc	Terminal capacitance Ct f=1MHz	NEP VR=70 V									
1,700 1101			λ=λρ	LSO 420 nm		CsI(TI) 540 nm	Isc 100 <i>lx</i>	Тур.	Max.	VR=70 V	VR=70 V	VR=70 V										
	(nm)	(nm)	(A/W)	(A/W)	(A/W)	(A/W)	(µA)	(nA)	(nA)	(times/°C)	(MHz)	(pF)	(W/Hz1/2)									
S3590-08		0 960	0.66	0.20	0.30	0.36	100	2	6	1.12	40	40	3.8 × 10 ⁻¹⁴									
S3590-09	340 to 1100		0.00	0.22	0.33	0.41	90						3.0 X 10 11									
S3590-18			0.65	0.28	0.34	0.38	100						7.6 × 10 ⁻¹⁴									
S3590-19														0.58	0.33	0.37	0.4	86	4	10		

PMT (H5783-04) specifications



According to manual at http://html.alldatasheet.com/html-pdf/62598/HAMAMATSU/H5783-04/407/1/H5783-04.html

Product Variations

Suffix Type No.	None	-01	-02	-03	-04	-06	-20	Output Type	Features
H5773	y e s	y e s	y e s	y e s	y e s	y e s	y e s	On-board	Low power consumption
H5783	y e s	y e s	y e s	y e s	y e s	y e s	y e s	Cable output	
H5773P	y e s	no	no	no	no	no	no	On-board	For photon counting
H5783P	y e s	no	no	no	no	no	no	Cable output	Low power consumption
H6779	y e s	y e s	y e s	y e s	y e s	y e s	y e s	On-board	Low ripple noise
H6780	y e s	y e s	y e s	y e s	y e s	y e s	y e s	Cable output	Fast settling time

Suffix	Spectral Response
None	300 nm to 650 nm
-01	300 nm to 850 nm
-02	300 nm to 880 nm
-03	185 nm to 650 nm
-04	185 nm to 850 nm
-06	185 nm to 650 nm
-20	300 nm to 900 nm

The suffix -06 type (synthetic silica window) has higher sensitivity than the -03 type below 300 nm in wavelength range.

Specifications

Parameter	H5773 / H5783 / H6779 / H6780 Series								
Suffix	None	-03, -06	-01, -04	-02	-20	_			
Input Voltage	+11.5 to +15.5								
Max. Input Voltage		V							
Max. Input Current	H5773 / H5783 Series: 9								
	H6779 / H6780 Series: 30								
Max. Output Signal Current	100								
Max. Control Voltage	+1.0 (Input impedance 100 k Ω)								
Recommended Control Voltage Adjustment Range +0.25 to +0.9									
Effective Area	♦ 8								